



Journey™ Series – Effortless

Wheelchair Cushion



Visco Elastic Memory Foam

Waterproof Breathable Layer




High Resiliency Pressure-Redistributing Foam

Waterproof Zipper

Waterfall Layer

Non-Skid Base



The coccyx, trochanters, ischial tuberosities, and the popliteal space are the areas that are at the greatest risk for skin breakdown and compromised circulation.

Product Specifications

HCPC E2601 (General Use)

Wound Risk Level: No History of Wounds or Risk Factors

Cushion Thickness / Weight Limit: 2" / 250lbs

Product Description

The Journey Effortless cushion is designed for individuals in need of a light-weight, general use cushion. Created with both the caregiver and user in mind, the Journey Effortless provides comfort for those with no risk of skin breakdown. A maintenance free solution, the Journey Effortless is the only Keen cushion with a nylon cover that can be easily wiped down for disinfecting/cleaning.

Models and Sizes Available (WIDTH x DEPTH x HEIGHT)

Model #	Size
F2 – 18 x 16	18" x 16" x 2" (6 cushions per case)

Comfort • Safety • Mobility

Toll-Free 1(866)330-5336

Keen

Copyright © 2007, Keen Mobility Company. All rights reserved. Rev 070312. All trademarks are property of their respective owner.

Seating 47