


Journey™ Series – Solid Seat Insert

Sling Seat Accessory



Are the hips to the back of the wheelchair? (continued from page 68)

If not, consider:

1. Evaluate for pelvic obliquity.
2. If a pelvic obliquity is present then carve a leg trough for the hip that lacks 90° hip flexion.
3. A flexible pelvic obliquity can be corrected using a solid seat insert or pelvic obliquity build-up.

JourneyTM Series – Solid Seat Insert

Sling Seat Accessory

Product Specifications

DO NOT PLACE INDIVIDUAL DIRECTLY ON INSERT

Weight Limit 600lbs

Product Description

The new Journey Solid Seat Insert is designed to offer a flat surface for wheelchairs with slinging or hammocking seats. Fully covered in non-skid material, this insert will stay in place. Use this insert with a Journey cushion to help eliminate internal rotation of the hips. Insert is properly labeled to help avoid putting an individual directly onto the insert. Insert must be used with a wheelchair cushion.

Models and Sizes Available (WIDTH x DEPTH)

Model #	Size
SSI – 16 x 16	16" x 16"
SSI – 16 x 18	16" x 18"
SSI – 18 x 16	18" x 16"
SSI – 18 x 18	18" x 18"
SSI – 20 x 18	20" x 18"

Comfort • Safety • Mobility

Toll-Free 1(866)330-5336

Keen[®]

Copyright © 2007, Keen Mobility Company. All rights reserved. Rev 070312. All trademarks are property of their respective owner.

Seating 71