

Liberty™ Series – Safety Mat

Fall Safety for the Bed



Easy for:

- Transport to other rooms
- Wheelchair-to-bed patient transfers
- Caregivers to store under beds



Tri-Fold Design

Non-Skid Bottom

Easy Grip Handle



- Tissue Load is the distribution of pressure, friction, and shear on the tissue.
- Positioning and support surfaces for individuals in bed are important to the management of tissue loads.
- Limit the amount of time the head of the bed is elevated since when the head of the bed is elevated, shear forces are exerted on the body.

Liberty™ Series – Safety Mat

Fall Safety for the Bed

Product Specifications

Purpose:	Fall Safety
Design:	Tri-Fold
Bottom:	Non-Skid & Nylon
Cushion Thickness:	2"

Product Description

The Liberty Safety Mat is designed to offer support for individuals in case of a fall in both bed and rehab environments. The Liberty Safety Mat is light weight, equipped with handles, tri-fold design for easy transport and an aesthetically pleasing dark green color. Designed with a non-skid bottom and a durable, tear-resistant nylon cover, the Safety Mat comes with 3 beveled edges. This allows for easy wheelchair and care provider access so that individual transfers are safe and simple.

Models and Sizes Available

Model #	Size (WIDTH x LENGTH x HEIGHT)
LSB2 – 64 x 36	64" x 36" x 2"

Comfort • Safety • Mobility

Toll-Free 1(866)330-5336

Keen

Copyright © 2007, Keen Mobility Company. All rights reserved. Rev 070312. All trademarks are property of their respective owner.

Bed Safety 27